

**BIG SAVINGS ON
SMALL INSTRUMENTS**

Valid till 31 December 2015

Hotplate Stirrers

Magnetic Stirrers

Vortex Mixers

Shakers

Thermoshakers

pH Meters

Pipettors



**Save BIG on Small
Laboratory Instruments!
See inside for details**

IKA[®]



VWR[®]
We Enable Science

Grant

SPRING 2015

HOTPLATE STIRRERS

SPECIAL OFFER



**Magnetic Stirrer with Hotplate, Ceramic Top
7 x 7" for up to 10L**

Technical Data

Stirring quantity max. (H2O) [l]	10
Speed range [rpm]	100 - 1500
Heating temperature range [°C]	50 - 500
Set-up plate material	ceramic
Set-up plate dimensions [mm]	180 x 180
Dimensions (W x H x D) [mm]	220 x 105 x 330
Protection class according to DIN EN 60529	IP 21

Cat. No.	Special Price
IKAA3581200	\$895.00

**DON'T MISS OUT ON THESE PRICES,
PLACE YOUR ORDER TODAY!**

**Best
Seller**



**Magnetic Stirrer with Hotplate, Aluminium
Plate for up to 20L**

Technical Data

Stirring quantity max. (H2O) [l]	20
Speed range [rpm]	0/50 - 1500
Heating temperature range [°C]	310
Set-up plate material	Aluminium alloy
Set-up plate dimensions [mm]	Ø 135
Dimensions (W x H x D) [mm]	160 x 85 x 270
Protection class according to DIN EN 60529	IP 42

Cat. No.	Special Price
IKAA20002620	\$1,290.00

IKA®

HOTPLATE STIRRERS

SUPER DEAL



**Analogue Magnetic Stirrer with Hotplate,
Aluminium Plate 7 x 7" for up to 2.5L**

Technical Data

Top Plate Dimensions 7 x 7" Units	7 x 7" (17.8 x 17.8cm)
Temperature Range Aluminium Top Plate Units	Ambient +5°C to 400°C
Speed Range	60 to 1600rpm
Maximum Capacity (H2O) 7 x 7" Units	2500mL
Overall Dimensions	14.8 x 9.9 x 4.3" (37.6 x 25.1 x 10.9cm)

Cat. No.

VWRI444-0602

Special Price

\$755.00



**Digital Magnetic Stirrer with Hotplate, 7 x 7"
Aluminium Plate for up to 2.5L**

Technical Data

Top Plate Dimensions 7 x 7" Units	7 x 7" (17.8 x 17.8cm)
Temperature Range Aluminium Top Plate Units	Ambient +5°C to 400°C
Speed Range	60 to 1600rpm
Timer	1 second to 160 hours
Maximum Capacity (H2O) 7 x 7" Units	2500mL
Overall Dimensions	14.8 x 9.9 x 4.3" (37.6 x 25 x 10.9cm)

Cat. No.

VWRI444-0632

Special Price

\$1,280.00

MAGNETIC STIRRERS

SPECIAL OFFER

**Best
Seller**



**Magnetic Stirrer with Digital Display
Colour Squid for up to 1L
Big Squid for up to 1.5L**

Technical Data - Colour Squid

Stirring quantity max. (H ₂ O) [l]	1
Speed range [1/min]	0 - 2500
Set-up plate dimensions [mm]	Ø 115
Dimensions (W x H x D) [mm]	145 x 45 x 160
Protection class according to DIN EN 60529	IP 54

Cat. No.

IKAA3671000

Special Price

\$270.00

Technical Data - Big Squid

Stirring quantity max. per stirring position (H ₂ O) [l]	1.5
Speed range [1/min]	0 - 2500
Set-up plate material	glass
Set-up plate dimensions [mm]	Ø 160
Dimensions (W x H x D) [mm]	180 x 40 x 195
Protection class according to DIN EN 60529	IP 54

Cat. No.

IKAA3672000

Special Price

\$365.00

**DON'T MISS OUT ON THESE PRICES,
PLACE YOUR ORDER TODAY!**



**Analogue Magnetic Stirrer with Ceramic Top,
7 x 7" for up to 10L**

Technical Data

Stirring quantity max. (H ₂ O) [l]	10
Speed range [1/min]	100 - 1500
Set-up plate material	ceramic
Set-up plate dimensions [mm]	200 x 200
Dimensions (W x H x D) [mm]	220 x 105 x 330
Protection class according to DIN EN 60529	IP 21

Cat. No.

IKAA3582400

Special Price

\$940.00

IKA®

MAGNETIC STIRRERS

**PRICE
REDUCED**



Analogue Magnetic Stirrer with Aluminium Plate, 7 x 7" for up to 2.5L

Technical Data

Top Plate Dimensions 7 x 7" Units	7 x 7" (17.8 x 17.8cm)
Temperature Range Aluminium Top Plate Units	Ambient +5°C to 400°C
Speed Range	60 - 1600 rpm
Maximum Capacity (H2O) 7 x 7" Units	2500mL
Overall Dimensions 7 x 7" Units	14.8 x 9.9 x 4.3" (37.6 x 25.1 x 10.9cm)

Cat. No.	Special Price
VWRI97042-632	\$470.00



Digital Magnetic Stirrer with Aluminium Plate, 7 x 7" for up to 2.5L

Technical Data

Top Plate Dimensions 7 x 7" Units	7 x 7" (17.8 x 17.8cm)
Temperature Range Aluminium Top Plate Units	Ambient +5°C to 400°C
Speed Range	60 to 1600rpm
Maximum Capacity (H2O) 7 x 7" Units	2500mL
Overall Dimensions 7 x 7" Units	14.8 x 9.9 x 4.3" (37.6 x 25 x 10.9cm)

Cat. No.	Special Price
VWRI97042-712	\$1,130.00

VORTEX MIXERS

SPECIAL OFFER



Analogue Vortex Mixer*

Technical Data	
Type of movement	orbital
Shaker diameter [mm]	4
Speed min (adjustable) [1/min]	500
Speed range [1/min]	0 - 2500
Operating mode	continuous operation
Touch function	yes
Working with microtiter plates	yes
Dimensions (W x H x D) [mm]	127 x 136 x 149

*Additional attachments available

Cat. No.	Special Price
IKAA3340000	\$495.00



Personal Vortex Mixer

Technical Data	
Speed	750 to 3000
Orbit	4
Maximum Tube Diameter	28.5
Capacity	1 up to 59 ml tube
Ambient Temperature Range	+ 4 to 40

Cat. No.	Special Price
GRANPV1	\$350.00

IKA®

SHAKERS



Multi Function 3D Rocker-Rotator

Technical Data	
Speed	1 to 100rpm
Tilt Angle	7°
Rocking Angle (Vibration Mode)	0 to 5° (step 1°)
Orbit Diameter	22mm
Timer for Orbital and Reciprocation Mode	0 to 250 sec
Timer for Vibration Mode	0 to 5 sec
Maximum Load	1kg
Display	2 x 16 character LCD

Cat. No.	Special Price
GRANPS-M3D	\$1,120.00

SUPER DEAL



Orbital Multi-Platform Shaker*

Technical Data	
Speed	20 to 50rpm
Orbit	20mm
Maximum Load	8kg
Timer, with Automatic Switch-Off	1 min to 96 hours
Motion Timer (Orbital/Reciprocal Modes)	0 to 250 sec
Motion Timer (Vibration Mode)	0 to 5 sec
Angle (Reciprocal Mode)	0 to 360°
Angle (Vibration Mode)	0 to 5°

* Platform sold separately

Cat. No.	Special Price
GRANPSU-20I	\$2,450.00

Grant

www.globalscience.co.nz

THERMOSHAKERS

**GREAT
DISCOUNT**



Heated Thermoshaker Ambient +5°C to 100°C*

Technical Data

Temperature Setting Range	+ 25 to 100
Temperature Control Range	+5°C above ambient to +100°C
Temperature Uniformity Over the Block	@ +37°C ± 0.1°C @ +60°C ± 0.2°C @ +100°C ± 0.2°C
Heat Up Speed RT to 100°C	4°C/min
Shaking Speed	250 to 1400
Orbit Diameter	2mm
Timer	1 min to 96 hours (1 min increment)

*Blocks sold separately

Cat. No.	Special Price
GRANPHMT	\$2,100.00



Heated and Cooled Thermoshaker Ambient -15°C to 100°C*

Technical Data

Temperature Setting Range	+ 4 to 100
Temperature Control Range	+15°C below ambient to +100°C
Temperature Uniformity Over the Block	@ +4°C ± 0.6°C @ +37°C ± 0.1°C @ +100°C ± 0.3°C
Heat Up Speed RT to 100°C	5°C/min (~ 15 mins) 1.8°C/min (~ 8 mins)
Shaking Speed	250 to 1400
Orbit Diameter	2mm
Timer	1 min to 96 hours (1 min increment)

*Blocks sold separately

Cat. No.	Special Price
GRANPCMT	\$2,200.00

Grant

pH METERS

**Best
Seller**



Seven Compact pH Meter & Probe Kit*

Great for high-end labs

Technical Data

Measuring Range	0.00...14.00 pH ± 0.1°C between 30°C and 130°C -2000...2000 mV 0°C...100°C
Resolution	0.001 / 0.01pH 0.1/1 mV / 0.1°C
Accuracy	± 0.002 pH / 0.2 mV / ± 0.1°C
Calibration	Max 5 points, 8 predefined to 1 user defined buffer groups
Input	BNC (Impedance > 3 x 10 ¹² Ω) Cinch / RCA (NTC30k Ω or PT1000 Reference: 2mm)
Temperature-compensation	ATC & MTC

*Kit includes meter, electrode, electrode arm & buffer sachets

Cat. No.	Special Price
METT30019029	\$1,885.00

GREAT PRICE



FiveEasy Compact pH Meter & Probe Kit*

For your basic laboratory pH needs

Technical Data

Measuring Range	0.00...14.00 pH ± 0.5°C between 20°C and 45°C -1999...1999 mV 0°C...100°C
Resolution	0.01 pH / 1mV / 0.1°C
Error Limits	± 0.01 pH / ± 1 mV / ± 0.5°C
Calibration	2 points, 4 predefined buffer groups
Input	BNC, (Impedance > 10e + 12Ω) Cinch / RCA (NTC 30 Ω) Reference: 2mm
Temperature-compensation	ATC & MTC

*Kit includes meter, electrode, electrode arm & buffer sachets

Cat. No.	Special Price
METT51302911	\$1,050.00

METTLER TOLEDO

ELECTRONIC PIPETTORS

**\$100 worth of FREE
Rainin tips with every
E4 Electronic Pipette**



Lite Touch System*

25% less force required

Cat. No.	Description	Unit	Price
RAIN17014484	Pipette 0.5-10µL E4 XLS Electronic with LTS Shaft	1 item	\$1,490.00
RAIN17014487	Pipette 2-20µL E4 XLS Electronic with LTS Shaft	1 item	\$1,490.00
RAIN17014483	Pipette 10-100µL E4 XLS Electronic with LTS Shaft	1 item	\$1,490.00
RAIN17014486	Pipette 20-200µL E4 XLS Electronic with LTS Shaft	1 item	\$1,490.00
RAIN17014488	Pipette 20-300µL E4 XLS Electronic with LTS Shaft	1 item	\$1,490.00
RAIN17014482	Pipette 100-1000µL E4 XLS Electronic with LTS Shaft	1 item	\$1,490.00
RAIN17014485	Pipette 200-2000µL E4 XLS Electronic with LTS Shaft	1 item	\$1,490.00
RAIN17012312	Pipette 500-5000µL E4 XLS Electronic with LTS Shaft	1 item	\$1,490.00
RAIN17012313	Pipette -1-10ML E4 XLS Electronic with LTS Shaft	1 item	\$1,490.00
RAIN17012314	Pipette 2-20ML E4 XLS Electronic with LTS Shaft	1 item	\$1,490.00

*LTS Tips Required

PLUS FREE TIPS!

Standard Shaft†

Cat. No.	Description	Unit	Price
RAIN17014491	Pipette 0.5-10µL E4 XLS Electronic with Standard Shaft	1 item	\$1,490.00
RAIN17014494	Pipette 2-20µL E4 XLS Electronic with Standard Shaft	1 item	\$1,490.00
RAIN17014490	Pipette 10-100µL E4 XLS Electronic with Standard Shaft	1 item	\$1,490.00
RAIN17014493	Pipette 20-200µL E4 XLS Electronic with Standard Shaft	1 item	\$1,490.00
RAIN17014495	Pipette 20-300µL E4 XLS Electronic with Standard Shaft	1 item	\$1,490.00
RAIN17014489	Pipette 100-1000µL E4 XLS Electronic with Standard Shaft	1 item	\$1,490.00
RAIN17014492	Pipette 200-2000µL E4 XLS Electronic with Standard Shaft	1 item	\$1,490.00
RAIN17012353	Pipette 500-5000µL E4 XLS Electronic with Standard Shaft	1 item	\$1,490.00
RAIN17012354	ELECTR. PIPETOR 1-10ML with Standard Shaft	1 item	\$1,490.00

†Fits most standard tip brands

PLUS FREE TIPS!

8 and 12 Multi channel LTS
and 6 and 8 Channel
Adjustable spacers LTS also
available, pricing on request.

RAININ

Pipetting 360°

RAININ MANUAL PIPETTORS



Lite Touch System*

25% less force required

Cat. No.	Description	Unit	Price
RAIN17014393	PIPETTE 0,1-2µL PIPET-LITE XLS LTS	1 item	\$500.00
RAIN17014388	PIPETTE 0,5-10µL PIPET-LITE XLS LTS	1 item	\$500.00
RAIN17014392	PIPETTE 2-20µL PIPET-LITE XLS LTS	1 item	\$440.00
RAIN17014384	PIPETTE 10-100µL PIPET-LITE XLS LTS	1 item	\$440.00
RAIN17014391	PIPETTE 20-200µL PIPET-LITE XLS LTS	1 item	\$440.00
RAIN17014405	PIPETTE 20-300µL PIPET-LITE XLS LTS	1 item	\$440.00
RAIN17014382	PIPETTE 100-1000µL PIPET-LITE XLS LTS	1 item	\$440.00
RAIN17014390	PIPETTE 200-2000µL PIPET-LITE XLS LTS	1 item	\$440.00
RAIN17011790	PIPETTE 500-5000µL PIPET-LITE XLS LTS	1 item	\$500.00
RAIN17011783	PIPETTE 1-10ML PIPET-LITE XLS LTS	1 item	\$500.00
RAIN17011788	PIPETTE 2-20ML PIPET-LITE XLS LTS	1 item	\$530.00

*LTS Tips Required

Standard Shaft†

Cat. No.	Description	Unit	Price
RAIN17014413	PIPETTE 0,1-2µL PIPET-LITE XLS with Standard Shaft	1 item	\$500.00
RAIN17014409	PIPETTE 0,5-10µL PIPET-LITE XLS with Standard Shaft	1 item	\$500.00
RAIN17014412	PIPETTE 2-20µL PIPET-LITE XLS with Standard Shaft	1 item	\$440.00
RAIN17014408	PIPETTE 10-100µL PIPET-LITE XLS with Standard Shaft	1 item	\$440.00
RAIN17014411	PIPETTE 20-200µL PIPET-LITE XLS with Standard Shaft	1 item	\$440.00
RAIN17014414	PIPETTE 20-300µL PIPET-LITE XLS with Standard Shaft	1 item	\$440.00
RAIN17014407	PIPETTE 100-1000µL PIPET-LITE XLS with Standard Shaft	1 item	\$440.00
RAIN17014410	PIPETTE 200-2000µL PIPET-LITE XLS with Standard Shaft	1 item	\$440.00
RAIN17011801	PIPETTE 500-5000µL PIPET-LITE XLS with Standard Shaft	1 item	\$500.00
RAIN17011795	PIPETTE 1-10ML PIPET-LITE XLS with Standard Shaft	1 item	\$500.00

†Fits most standard tip brands



Cat. No.	Description	Unit	Price
RAIN17014406	KIT STARTER PIPET LITE XLS	1 KIT	\$1,260.00
RAIN17014415	STARTER KIT PIPETTE STANDARD SHAFT	1 KIT	\$1,260.00

Starter Kit Includes 20XLS, 200XLS, 1000XLS, 1 box of Tips for each and hang ups for each.

VWR COLLECTION PIPETTORS

3 YEAR WARRANTY



Ergonomically Designed, High Accuracy Precision Pipettes

These new VWR Standard Line pipettes provide accuracy and precision you can rely on. The ergonomic features deliver comfortable pipetting even during extensive periods of use.

- Lightweight
- Soft spring system for smooth, effortless pipetting
- Easily adjustable ejector height system to accommodate virtually all brands of tips
- UV resistant and fully autoclavable
- Easy and quick recalibration system
- Colour coded shaft for easy

- volume identification
- 5 and 10 ml shaft protected by filter

Multichannel Standard Line pipettes feature a revolutionary suspension system which allow the shafts to retract slightly when pressed against a row of pipette tips. Each shaft moves independently, resulting in an even sealing force being applied to all 8 or 12 tips. The suspension system not only ensures secure tip sealing but also eliminates tips becoming jammed. The 360° manifold on multichannels guarantees optimal working comfort.

Cat. No	VOLUME RANGE (µl)	Accuracy (%)	Imprecision (%)	Price
Single channel				
VWRI613-2891	0,2 - 2	±12,0 - ±1,5	±6,0 - ±0,7	\$196.00*
VWRI613-2892	0,5 - 10	±4,0 - ±0,5	±2,8 - ±0,4	\$196.00*
VWRI613-2893	2 - 20	±3,0 - ±0,8	±1,5 - ±0,3	\$196.00*
VWRI613-2894	10 - 100	±1,6 - ±0,8	±0,8 - ±0,2	\$196.00*
VWRI613-2895	20 - 200	±1,2 - ±0,6	±0,6 - ±0,2	\$196.00*
VWRI613-2896	100 - 1000	±1,6 - ±0,6	±0,4 - ±0,15	\$196.00*
VWRI613-2897	500 - 5000	±1,2 - ±0,5	±0,5 - ±0,15	\$220.00*
VWRI613-2898	1000 - 10000	±2,5 - ±0,5	±0,6 - ±0,2	\$220.00*
8 channels				
VWRI613-2899	0,5 - 10	±10 - ±2,0	±8,0 - ±1,2	\$490.00*
VWRI613-2900	5 - 50	±4,0 - ±1,6	±2,5 - ±0,6	\$490.00*
VWRI613-2901	20 - 200	±3,0 - ±1,0	±1,5 - ±0,6	\$490.00*
VWRI613-2902	50 - 300	±1,6 - ±1,0	±1,5 - ±0,6	\$490.00*
12 channels				
VWRI613-2903	0,5 - 10	±10 - ±2,0	±8,0 - ±1,2	\$595.00*
VWRI613-2904	5 - 50	±4,0 - ±1,6	±2,5 - ±0,6	\$595.00*
VWRI613-2905	20 - 200	±3,0 - ±1,0	±1,5 - ±0,6	\$595.00*
VWRI613-2910	50 - 300	±1,6 - ±1,0	±1,5 - ±0,6	\$595.00*

*Contact us for a further discount on a bulk purchase of 3 or more pipettes

IKA®

VWR 
We Enable Science

RAININ

Pipetting 360°

METTLER TOLEDO

Grant

www.globalscience.co.nz

Prices Valid Till 31/12/2015 | Prices exclude GST & Freight

Global Science a VWR Company

241 Bush Road, Albany 0632,
Auckland, New Zealand
Tel: 64 9 443 5867
Fax: 64 9 444 7314
E-mail: sales@globalscience.co.nz
Website: www.globalscience.co.nz

PRICES VALID TILL 31/12/2015.
SEE OUR WEBSITE FOR UP TO DATE
PRICING AND STOCK AVAILABILITY.
PRICES EXCLUDE GST AND FREIGHT
