

PALL LIFE SCIENCES SENTINO FILTER DISPENSER

Simple Design with
reliable dispense



NEW Sentino™ Filter Dispenser

Easy to load — Simple loading and unloading of refill pack at any point during use.

No jamming — Guide plates assist loading and keep refill web in alignment during use to avoid jamming of the instrument.

Soft touch keypad — Aseptically dispense a disc of membrane at the touch of a button.

Designed to hold 200 membrane discs — Saves time by eliminating interruptions to reload dispenser.

Compact design — Small footprint frees valuable bench top space and provides flexibility in arranging workspace for optimal efficiency and workflow.

Miniaturised refill pack — Sized to minimize storage space resulting in less packaging material waste.

Portable — Use plugged into main power source or charge battery for use anywhere. Indicator lights let you know when the battery charge is running low.

Easy to clean — Sanitize with standard laboratory disinfectants.

Ordering Information

Part Number	Pall PN	Description	Pkg
PALLCA30617-196	13182	Sentino Filter Dispenser	EA
Select Your Membrane Disk Type			
PALLCA30617-200	68123	0.2 µm, Supor® membrane	1000/pkg
PALLCA30617-198	68121	0.45 µm, GN-6 Metrice® membrane	1000/pkg
PALLCA30617-202	68124	0.45 µm, Metrice Black membrane	1000/pkg
PALLCA30617-204	68125	0.8 µm, Metrice Black membrane	1000/pkg