



INTRODUCING THE ALL NEW BUCHI ROTARY EVAPORATION FAMILY

R-300 System

**Your fully automated distillation
evaporation system**



- Centrally controlled vacuum pump
- Option of large (5L) or small (1L) bath
- Foam detection
- Automatic Distillation
- Remote monitoring via the Büchi phone app or wireless system
- Featuring the Büchi simple flask snap lock system
- Option of the I-300 or I-300 Pro to optimise your distillation process
- Choice of auto or manual lift function

CONTACT US FOR A QUOTE TODAY

Choose Your Controller



R-300 - your fully automated distillation evaporation system

R-300 Interface Controllers - Central Control of System Components

I-300



- Remote monitoring with Buchi phone app including push notifications and live view
- Multiple mounting position options
- Unattended operation with automatic distillation



I-300 Pro



- Central 7" touch screen control panel
- Guided operation with programmable methods
- Remote monitoring with Buchi phone app including push notifications and live view
- Multiple mounting position options



R-100 System

Your Economical Evaporation Solution



- Reliable operation due to the robust design
- Smooth manual lift mechanism for quick ergonomic operation
- Easy flask handling using the multi-functional combi-clip
- Easy operation by integrating an interface
- Flexibility to add vacuum pump and chiller

Choose Your Controller


R-100 Interface Controller - Functional, Easy to Use Design

I-100



- Convenient, easy to use design
- Robust and durable
- User friendly self-explanatory interface
- Highly reliable due to centrally controlled process parameters

Newly Designed Vacuum Source V-300



- Near silent operation
- Flow rate of 1.8m³/h
- Ultimate vacuum of 5mbar
- Optimal vacuum pump for the R-300 Rotovapor
- Unattended operation in combination with the I-300 or I-300pro and optional level sensor
- High chemical resistance
- Easy visual detection of condensate and particles through the transparent front panel

CONTACT US FOR A QUOTE TODAY

PRICED FROM \$3,900.00



Economical Vacuum V-100



- Small footprint to save space
- Multiple application use including evaporation, drying and concentrating and suction straining and filtration
- Light weight compact design
- Low maintenance system
- High chemical resistance

PRICED FROM \$2,400.00

Additional Accessories CHILLERS

F-305, F-308, F-314 F-100



**GO TO WWW.GLOBALSCIENCE.CO.NZ
FOR THE LATEST SPECIAL OFFERS**