

METTLER TOLEDO TOUCH SCREEN BALANCES

Introducing the new Touch Screen Balance Range from Mettler Toledo

With built in glove mode for touch screen operation
(operable through cotton, silicone and latex)

METTLER TOLEDO

MS-TS Range Features

- MS-TS range features a 7" colour touch screen versatile modular format
- Robust design with full metal housing
- Straight forward intuitive interface
- Visual level bubble on the front of the balance
- On screen instructions guide to level the balance
- Fully automatic time and temperature controlled internal calibration
- Min weigh functions for easy and reliable results

**DAILY TASKS BECOME EASIER
WITH THE OPTIMISED DISPLAY**

**EASY SCROLL FUNCTION
JUST LIKE A SMART PHONE**



Options include:
4 decimal place options: 100g- 300g
3 decimal place options: 200g-1kg
2 decimal place options: 800g-12kg
1 decimal place options: 3kg-32kg

**QuickLock draft
shield dismantles
in seconds
without any tools
and all parts are
dishwasher proof.**

**Safe, easy
cleaning with
smooth rounded
surfaces.**

METTLER TOLEDO TOUCH SCREEN BALANCES

Big Features Go Compact

With built in glove mode for touch screen operation
(operable through cotton, silicone and latex)

METTLER TOLEDO

ML-T Range Features

- Easy portable design with robust bottom metal housing
- Small space saving footprint
- Visual level bubble on the front of the balance
- On screen instructions guide to level the balance
- Fully automatic time and temperature controlled internal calibration
- Min weigh functions for easy and reliable results



**Featuring a 4.3"
colour touch
screen.**

**Optional battery
operation.**

Options include:

4 decimal place options: 100g- 300g

3 decimal place options: 200g-1kg

2 decimal place options: 800g-12kg

1 decimal place options: 3kg-32kg

**FOR THE FULL RANGE OF
WHAT WE OFFER AND A
QUOTE CALL US TODAY**